



Standard Fabrics:



Antler



Elk



Lakefront



Chestnut

This contemporary plush fabric comes in beautiful, rich colors. It features our exclusive SofTouch™ with Stain Defense coating.

Crossover fabrics are unavailable.

Eirene
PR202



Comfort Zone™
by GOLDEN

Inspired by the moon, the DeLuna™ Series of luxury power lift recliners is designed to offer you the most calming comfort and restorative sleep. The Eirene is medium-sized with a back style that provides maximum comfort in a tasteful button design. It features a supportive seating pad, an easy grip armrest design, and a USB port on the hand control for your phone or tablet. Available in luxurious new fabrics to match any décor.

Standard Features:

- Easy Grip Armrest
- Headrest Cover
- USB Charging Port



Quick Ship?

Check real-time inventory at webportalgt.com. (Access available for Golden retailers only.)

DeLUNA Series PR202

Model Number	PR202-MED
Size	Medium
Weight Capacity	375 LBS
Positioning	3
Comfort Zone(s)	1
Twilight Technology	NO
Chaise Seat	NO
Pocketed Coil Seat	NO
Back Type	Button
Removable Back Type	KD
Head Rest Cover	YES
Arm Rest Covers	NO
Overall Width	32"
Overall Height	43"
Width Between Arms	21"
Floor to Top of Seat	20.5"
Seat Depth	19.5"
Seat to Top of Back	28"
Seat to Top of Arm	6.5"
Distance Required from Wall Reclined	20"
Extension Reclined from Back of Seat to End of Foot Rest	43"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES
Warranty on Electrical Parts	3 YR
Pro Rated Years 4 Thru 7	YES
Warranty on Mechanical Labor	3 YR
Warranty on Electrical Labor	1 YR
Hand Control Routing	Inside Arm
Extra Pocket Option	NO
Heat & Massage Option	NO
Articulating Headrest Option	NO
Adjustable Lumbar Option	NO
AutoFlex Option	NO
Foot Rest Extension Option	NO
Straight Lift Option	YES
USB Charging Port	Hand Control
Battery Backup	YES
Place of Origin	China



Shown in
Lakefront

SmartTek Diagnostic System STANDARD on all Golden models.

Designed and engineered in Old Forge, PA. Assembled in China using Golden components.

Comfort Zone™ by GOLDEN



Comfort Zone 1 provides ergonomic placement of your lower body allowing your muscles to relax for long-lasting comfort. Comfort Zone 1 can be adjusted using the two buttons on the hand control.



WE PATENTED COMFORT™



Scan with your smartphone camera to learn more about Comfort Zone by Golden.



Phone: 800-624-6374
Fax: 800-628-5165
www.goldentech.com

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Complies with
UL 73 and
CSA C22.2 No. 68



www.facebook.com/goldentechnologies

SHT-PR202-042221