

Fly-Lite Aluminum Transport Chair

HCPCS: E1038



- Weighs only 16.2 lbs without footrests 18.9 lbs with footrests
- Aluminum frame is lightweight and strong
- Comes with carry pocket on backrest
- Comes with seat belt for added safety
- Attractive plaid style nylon upholstery (Figure A)
- Comes with handy carry pocket on back (Figure B)
- Tool-free, height-adjustable swing-away footrests come standard
- Optional leg rests available (Model# STDELR-TF)
- 8" rear wheel 6" casters in front
- Back folds down for storage and transport and features deluxe back release (Figure C)
- Comes with aluminum rear wheel locks
- Convenient cup holder can easily attach to chair Sold separately (Model# STDS1040S) (Figure D)
- Weight Capacity: 300 lbs

| PRODUCTS | Description |
|-----------|------------------------------------|
| DFL19-BL | Fly-LiteTransportChair,19"Bl,1/cs |
| DFL19-BLK | Fly-LiteTransportChair,19"Blk,1/cs |
| DFL19-RD | Fly-LiteTransportChair,19 Red,1/cs |

| SPECIFICATIONS | Description |
|-------------------------------|--------------------------|
| Back Height | 18" |
| Carton Shipping Weight | 23.8 lbs |
| Overall Height | 36" |
| Overall Length w/ Riggings | 33" |
| Overall Width | 22" (open) 9" (closed) |
| Seat Dimensions | 19" (W) x 15.5" (D) |
| Seat To Floor Height | 18" |
| Weight Capacity | 300 lbs |
| Weight With Front Riggings | 18.9 lbs |
| Weight Without Front Riggings | 16.2 lbs |

| ACCESSORIES | Description |
|-------------|---------------------------------|
| STDS1040S | Cup Holder, Swivel,1BX,6/CS, BX |

| PARTS | Description |
|------------|----------------------------------|
| STDS2A1807 | SEAT RAIL PAD,DFL,NITRO,1/EA, EA |
| STDS2A1803 | Cross Frame Assist DFL,1/EA, EA |

| PARTS | Description |
|------------|---|
| STDS2A184R | BRAKE ASSY RIGHT,DFL,1/EA, EA |
| STDS2A184L | BRAKE ASSY LEFT,DFL,1/EA, EA |
| STDS2A1810 | FRONT WHEEL,DFL,TS19,1/ea, EA |
| STDS2A1806 | FRONT AXLE BOLT,DFL,1 PR, PR |
| STDS2A1811 | REAR AXLE BOLT,DFL,1 PR, PR |
| STDS2A1801 | REAR HUB CAP,DFL,1 PR, PR |
| STDS2A1800 | Rear Wheel, DFL, TS19,1/ea, EA |
| TRB-1 | BackSupportReleaseLever./Trans,1/ea, EA |
| STDS2A181L | HANDLE BAR,LEFT,BLACK,DFL,1/EA, EA |
| STDS2A182L | HANDLE BAR,LEFT,BLUE,DFL,1/EA, EA |
| STDS2A1891 | Back Upholstery, Black/DFL,1/ea, EA |
| STDS2A1892 | Back Upholstery, SilverPlaid,DFL,1/EA, EA |
| STDS2A1893 | BACK UPHOLSTERY Blue,DFL,1/EA, EA |
| STDS001014 | Folding Cane Hindge DFL,1/ea, EA |
| STDS2A1824 | Arm Pad,Black,1/ea, EA |
| STDS2A1804 | Arm Pad,Black Plaid,DFL,1ea, EA |
| STDS2A1814 | Arm Pad,Blue Plaid,DFL,1ea, EA |
| DFLSFBK | DFL Footrests, Black, 1pr/cs, CS |