

NasalFit Deluxe EZ CPAP Mask

HCPCS: A7032 A7034 A7035-1,A7032,A7034,A7035



FIGURE A



FIGURE B



FIGURE C

- The NasalFit Deluxe EZ CPAP Mask has a soft silicon cushion, comfortable forehead plate with silicon pad, and a 360° rotational elbow, which are designed to increase comfort, compliance, and mask stability
- Adjustable headgear, easy to assemble cushions, and a lightweight streamlined design help to ensure proper fit and minimal leak during sleep
- By providing a mask option without headgear, a more cost effective option is available based upon the reimbursement schedule
- Latex Free

PRODUCTS	Description
100NDEL	NasalFit Deluxe EZ CPAP Mask, LG, 1/ea
100NDEM	NasalFit Deluxe EZ CPAP Mask, MD, 1/ea
100NDES	NasalFit Delux EZ CPAP Mask, SM, 1/ea
ACCESSORIES	Description
100NDEH	Replacement Headgear for 100NDE, 1/ea, EA
100NDES-CUSHION	Replacement Cushion EZ Nasal, SM, 1/ea, EA
100NDEM-CUSHION	Replacement Cushion EZ Nasal, MD, 1/ea, EA
100NDEL-CUSHION	Replacement Cushion EZ Nasal, LG, 1/ea, EA
100NDEM-NH	NasalFit Deluxe EZ CPAPMask,NoHDG,MD,1/ea, EA
100NDEL-NH	NasalFit Deluxe EZ CPAPMask,NoHDG,LG,1/ea, EA
100NDES-NH	NasalFit Delux EZ CPAPMask,NoHDG, SM,1/ea, EA